


Makromed

TECHNOLOGIES



**EXTRACTION
FORCEPS
(ENGLISH
PATTERN)**

SURGICAL INSTRUMENTS
CATALOGUE



APPLICATION CHART FOR

EXTRACTION FORCEPS





Fig	Application
1	Upper centrals, canines and roots
2	Upper laterals
7	Upper premolars and roots
13	Lower premolars
17	Upper molars (Right)
18	Upper molars (Left)
19	Upper wisdoms
22	Lower molars (Hawk's bill)
29	Upper roots
29S	Upper roots
30	Upper roots
37	Children's upper incisors and canines
39	Children's upper molars
44N	Upper roots
51	Upper wisdoms (Bayonet)
67	Upper wisdoms
73	Lower molars (Hawk's bill)
73S	Lower molars (Hawk's bill)
74	Lower roots
74N	Small lower roots and anteriors
75	Lower premolars
76	Upper roots (Read)
76S	Upper roots (Read)
76N	Small upper roots (Read)

Fig	Application
79	Lower wisdoms
86	Lower molars (Cowhorn)
87	Lower molars (Cowhorn)
89	Upper molars (Right)(Upper Cowhorns)
90	Upper molars (Left)(Upper Cowhorns)
94	Upper molars (Right)
95	Upper molars (Left)
101	Upper molars and roots
107	Upper canines
113	Upper roots (Martin)
123	Children's lower anteriors, premolars, roots
136	Upper premolars and roots (Guy's)
137	Upper anteriors and roots (Guy's)
138	Children's upper anteriors, premolars, roots
147	Upper roots (Gibb)
157	Children's upper molars
158	Children's upper molars
159	Children's upper premolars
160	Children's lower molars (Hawk's bill)
161	Children's lower molars (Hawk's bill)
162	Children's lower anteriors and roots
163	Children's upper anteriors
3	Rongeurs forceps

HANDLE OPTIONS

SMALL MODEL



LARGE MODEL



FINGER-SHAPED



DIAMOND TIPPED





450-01 Fig. 01
Upper centrals and canines, wide



450-02 Fig. 02
Upper laterals and canines



450-03 Fig. 03
Upper crowded incisors & canines



450-04 Fig. 04
Lower incisors and canines



450-07 Fig. 07
Upper Premolars



450-08 Fig. 08
Lower premolars either side



450-13 Fig. 13
Lower premolars either side



450-13S Fig. 13S
Lower premolars, for children




450-17 Fig. 17
Upper molars, right



450- 18 Fig. 18
Upper molars, left



450- 18A Fig. 18A
Upper molars either side



450- 18B Fig. 18B
Upper molars



450- 19 Fig. 19
Upper third molars & upper
wisdoms either side



450- 20 Fig. 20
Lower third molars



450- 21 Fig. 21
Lower molars



450- 22 Fig. 22
Lower molars



450- 22S Fig. 22S
Lower molars, for children



450- 23 Fig. 23
Lower molars, right



450- 24 Fig. 24
Lower molars, left



450- 29 Fig. 29
Upper roots and incisors



450- 29N Fig. 29N
Root of Upper Incisors



450- 29S Fig. 29S
Upper roots, for children



450- 30 Fig. 30
Upper roots



450- 31 Fig. 31
Lower roots & incisors with
Serrated tips




450- 32 Fig. 32
Lower molars




450- 33 Fig. 33
Lower roots



450- 33A Fig. 33A
Lower roots



450- 33C Fig. 33C
Lower roots



450- 33L Fig. 33L
Lower roots



450- 33M Fig. 33M
Lower roots with long beak



450- 34 Fig. 34
Incisors and Bicuspid, deep gripping



450- 34N Fig. 34N
Incisors and Bicuspid, deep gripping



450- 35 Fig. 35
Incisors and Bicuspid, deep gripping



450- 35N Fig. 35N
Incisors and Bicuspid, deep gripping



450- 36 Fig. 36
Incisors and Bicuspid, deep gripping



450- 37 Fig. 37
Upper incisors, for children



450- 38 Fig. 38
Lower incisors and
canines, for children



450- 39 Fig. 39
Upper molars, for children



450- 39A Fig. 39A
Upper molars, for children



450- 39L Fig. 39L
Upper molars left, for children



450- 39R Fig. 39R
Upper molars right, for children



450- 40 Fig. 40
Lower molars, for children



450- 41 Fig. 41
Upper roots



450- 44 Fig. 44
Upper roots



450- 44F Fig. 44F
Upper roots



450- 45 Fig. 45
Lower roots



450- 45F Fig. 45F
Lower roots



450- 46 Fig. 46
Lower roots



450- 46L Fig. 46L
Lower roots



450- 46LX Fig. 46LX
for very fine roots



450- 46N Fig. 46N
Bicuspid, deep gripping




450- 49 Fig. 49
Upper roots



450- 50 Fig. 50
Upper roots




450- 51 Fig. 51
Upper roots



450- 51A Fig. 51A
Upper roots



450- 51C Fig. 51C
Upper roots



450- 51L Fig. 51L
Upper roots



450- 51M Fig. 51M
Especially long beaks



450- 52 Fig. 52
Upper roots



450- 52A Fig. 52A
Upper teeth & Resection



450- 53 Fig. 53
Upper roots




450- 54 Fig. 54
For separating upper molars




450- 55 Fig. 55
For separating upper molars



450- 56 Fig. 56
For separating lower molars



450- 65L Fig. 65L
Upper molars, left



450- 65R Fig. 65R
Upper molars, right




450- 67 Fig. 67
Upper third molars and
upper wisdoms either side



450- 67A Fig. 67A
Upper third molars and upper
wisdoms either side



450- 67LX Fig. 67LX
Upper third molar, left



450- 67RX Fig. 67RX
Upper third molar, right



450- 67N Fig. 67N
Third molar, deep-gripping



450- 68 Fig. 68
Lower roots



450-73 Fig. 73
Lower molars



450-73L Fig. 73L
Lower molars



450-73S Fig. 73S
Lower molars



450-74 Fig. 74
Lower roots and incisors



450-74D Fig. 74D
Lower roots



450-74M Fig. 74M
Lower roots and incisors



450-74N Fig. 74N
Lower roots and incisors




450-75 Fig. 75
Lower incisors and premolars



450-75A Fig. 75A
Lower incisors and premolars




450- 76 Fig. 76
Upper roots



450- 76N Fig. 76N
Upper roots



450- 76S Fig. 76S
Upper roots



450- 77 Fig. 77
Upper roots



450- 79 Fig. 79
Lower third molars and
lower wisdoms either side



450- 79A Fig. 79A
Lower third molars



450- 79N Fig. 79N
Lower Wisdoms



450- 86 Fig. 86
Lower molars decayed
or broken down crowns



450- 86B Fig. 86B
Lower molars



450- 86C Fig. 86C
Lower molars decayed or
broken down crowns



450- 87 Fig. 87
Lower molars decayed or
broken down crowns




450- 89 Fig. 89
Upper molars, right




450- 90 Fig. 90
Upper molars, left




450- 91 Fig. 91
For separating lower molars



450- 92 Fig. 92
Upper molars, right



450- 93 Fig. 93
Upper molars, left



450- 94 Fig. 94
Upper molars, right




450- 95 Fig. 95
Upper molars, left



450- 96 Fig. 96
Upper molars, right




450- 97 Fig. 97
Upper molars, left




450- 99 1/4 Fig. 99 1/4
Lower molars



450- 99 Fig. 99
Lower molars




450- 101 Fig. 101
Upper roots




450- 101A Fig. 101A
Upper roots




450- 107 Fig. 107
Upper Canines



450- 108 Fig. 108
Upper roots



450- 110 Fig. 110
Upper roots



450- 111 Fig. 111
Upper roots



450- 112 Fig. 112
Upper roots



450- 113 Fig. 113
Upper roots



450- 114 Fig. 114
Lower roots & centrals



450- 121 Fig. 121
Upper wisdoms



450- 122 Fig. 122
Lower centrals & roots



450- 123 Fig. 123
Lower teeth & roots,
for children




450- 130 Fig. 130
Upper wisdoms 190 mm



450- 133 Fig. 133
for very fine roots



450- 136 Fig. 136
Upper Centrals,
Premolars & roots




450- 136H Fig. 136H
Upper roots



450- 137 Fig. 137
Lower centrals & roots



450- 138 Fig. 138
Upper teeth & roots,
for children



450- 139 Fig. 139
Upper roots



450- 143 Fig. 143
Lower roots & centrals



450- 144 Fig. 144
Lower roots & centrals



450- 145 Fig. 145
Lower roots



450- 146 Fig. 146
Root fragments,
deep gripping




450- 147 Fig. 147
Upper roots



450- 148 Fig. 148
Lower molars



450- 151 Fig. 151
for very fine roots



450- 153 Fig. 153
Root fragments,
deep gripping



450- 155 Fig. 155
Lower roots



450- 157 Fig. 157
Upper molars,
for children



450- 158 Fig. 158
Upper premolars & molars,
for children



450- 159 Fig. 159
Upper premolars,
for children



450- 160 Fig. 160
Lower molars,
for children



450- 161 Fig. 161
lower molars, for children



450- 162 Fig. 162
Lower roots & centrals,
for children



450- 163 Fig. 163
Upper centrals, for children



450- 166 Fig. 166
Upper universal application



450- 167 Fig. 167
Lower teeth, universal pattern



450- 168 Fig. 168
Upper central & lateral teeth




450- 169 Fig. 169
Premolars, lower canines
& all roots



450- 171 Fig. 171
Wisdoms & molar roots



450- 173 Fig. 1
Lingually inclined roots and
wisdom teeth



450- 180 Fig. 180
Special for upper centrals and
canines deep gripping



450- 181 Fig. 181
Special for upper premolars,
deep gripping




450- 182 Fig. 182
Lower centrals canines and
premolars deep gripping



450- 185 Fig. 185
Lower roots



450- 187 Fig. 187
Lower roots fragments



450- 189 Fig. 189
Fragments upper roots



450- 190 Fig. 190
Fragments upper roots



450- 221 Fig. 221
Universal upper roots 150 mm



450- 223 Fig. 223
Universal for lowr roots



450- 350 Fig. 350
Lower molars & wisdoms left



450- 351 Fig. 351
Lower molars & wisdoms right

Makromed
TECHNOLOGIES



✉ info@makromed.tech | sales@makromed.tech

🌐 www.makromed.tech

