



## MOUTH GAGS & RETRACTORS

Mouth Gags  
gags bouche dentaires  
gagurile gura  
gag bocca  
gags boca

**Molt Mouth Gag (Adult)**  
with Silicone Tubing

| Catalog # | Size  |
|-----------|-------|
| -670-01   | 14 cm |



**Molt Mouth Gag (Child)**  
with Silicone Tubing

| Catalog # | Size    |
|-----------|---------|
| MM-670-02 | 11.4 cm |



**Molt Mouth Gag (Pediatric)**  
with Silicone Tubing

| Catalog # | Size    |
|-----------|---------|
| MM-670-03 | 10.7 cm |



**Roser Koenig Mouth Gag**

| Catalog # | Size   |
|-----------|--------|
| MM-670-04 | 160 mm |
| MM-670-05 | 185 mm |



**Denhardt Mouth Gag**  
with Grip Lock

| Catalog # | Size   |
|-----------|--------|
| MM-670-06 | 155 mm |



Mouth Gags  
 gagurile gura  
 bouche prop  
 bocca bavaglio  
 boca gag

### Fergusson-Ackland Mouth Gag

| Catalog # | Size   |
|-----------|--------|
| MM-670-07 | 145 mm |



### Heister Mouth Gag

| Catalog # | Size   |
|-----------|--------|
| MM-670-08 | 130 mm |



### Jennings Mouth Gag

| Catalog # | Size   |
|-----------|--------|
| MM-670-09 | 110 mm |
| MM-670-10 | 130 mm |
| MM-670-11 | 150 mm |



**Weitlaner Self Retaining Retractor**  
2 x 3 Sharp Prongs

| Catalog # | Size         |
|-----------|--------------|
| MM-990-01 | 11 cm (4 ¼") |




**Weitlaner Self Retaining Retractor**  
3 x 4 Sharp Prongs

| Catalog # | Size   |
|-----------|--------|
| MM-990-02 | 155 mm |



**Tobold Tongue Depressor**

| Catalog # | Size                   |
|-----------|------------------------|
| MM-990-03 | 15 cm (6"), 88 x 25 mm |



**Andrew Tongue Depressor**

| Catalog # | Size           |
|-----------|----------------|
| MM-990-04 | 14.5 cm (5 ¾") |



**Kroenig Tongue Depressor**

| Catalog # | Size       |
|-----------|------------|
| MM-990-05 | 15 cm (6") |



**Standard Tongue Depressor**

| Catalog # | Size           |
|-----------|----------------|
| MM-990-06 | 14.5 cm (5 ¾") |



Retractors  
retractores  
Ecarteurs  
retractoare  
Divaricatori

### Hartmann Tongue Depressor

| Catalog # | Size   |
|-----------|--------|
| MM-990-07 | 140 mm |



### Bruenings Tongue Depressor

| Catalog # | Size                        |
|-----------|-----------------------------|
| MM-990-08 | 19 cm (7 ½"), 19 mm - 15 mm |



### Minnesota Cheek Retractor

| Catalog # | Size           |
|-----------|----------------|
| MM-990-09 | 14.5 cm (5 ¾") |



### Roux Retractor

Set of Fig. 1, Fig. 2 and Fig. 3

| Catalog # | Size |
|-----------|------|
| MM-990-10 | N/A  |



### Sternberg Retractor

| Catalog # | Size           |
|-----------|----------------|
| MM-990-11 | 12 cm (4 ¾")   |
| MM-990-12 | 13 cm (5")     |
| MM-990-13 | 13.5 cm (5 ¼") |






Retractors  
retractores  
Ecarteurs  
retractoare  
Divaricatori

### Langenbeck Retractor


| Catalog # | Size, Blade Size           |
|-----------|----------------------------|
| MM-990-14 | 22.5 cm (8 ¾"), 30 x 11 mm |
| MM-990-15 | 22.5 cm (8 ¾"), 30 x 14 mm |
| MM-990-16 | 22.5 cm (8 ¾"), 30 x 16 mm |
| MM-990-17 | 22.5 cm (8 ¾"), 40 x 11 mm |
| MM-990-18 | 22.5 cm (8 ¾"), 50 x 11 mm |



### Langenbeck Retractor

#### Hollow Handle

| Catalog # | Size, Blade Size           |
|-----------|----------------------------|
| MM-990-19 | 22.5 cm (8 ¾"), 30 x 11 mm |
| MM-990-20 | 22.5 cm (8 ¾"), 30 x 14 mm |
| MM-990-21 | 22.5 cm (8 ¾"), 30 x 16 mm |
| MM-990-22 | 22.5 cm (8 ¾"), 40 x 11 mm |
| MM-990-23 | 22.5 cm (8 ¾"), 50 x 11 mm |



### Kocher Retractor

| Catalog # | Size, Blade Size       |
|-----------|------------------------|
| MM-990-24 | 23 cm (9"), 41 x 18 mm |
| MM-990-25 | 23 cm (9"), 61 x 20 mm |
| MM-990-26 | 23 cm (9"), 61 x 25 mm |
| MM-990-27 | 23 cm (9"), 75 x 30 mm |
| MM-990-28 | 23 cm (9"), 75 x 40 mm |




### Buchs Retractor

| Catalog # | Size, Blade Size         |
|-----------|--------------------------|
| MM-990-29 | 20.5 cm (8"), 38 x 18 mm |



### Middeldorf Retractor

| Catalog # | Size, Blade Size           |
|-----------|----------------------------|
| MM-990-30 | 21.5 cm (8 ½"), 15 x 17 mm |
| MM-990-31 | 23 cm (9"), 22 x 22 mm     |
| MM-990-32 | 24.5 cm (9 ¾"), 26 x 31 mm |



Retractors  
retractores  
Ecarteurs  
retractoare  
Divaricatori

### Halsted-Mosquito Heamostatic Forcep Straight

| Catalog # | Size           |
|-----------|----------------|
| MM-880-01 | 12.5 cm (5")   |
| MM-880-02 | 14.5 cm (4 ¾") |



### Halsted-Mosquito Heamostatic Forcep Curved

| Catalog # | Size         |
|-----------|--------------|
| MM-880-03 | 12 cm (4 ¾") |
| MM-880-04 | 14 cm (5 ½") |



### Halsted-Mosquito Heamostatic Forcep Straight, 1 x 2 Teeth

| Catalog # | Size           |
|-----------|----------------|
| MM-880-05 | 12.5 cm (5")   |
| MM-880-06 | 14.5 cm (5 ¾") |



### Halsted-Mosquito Heamostatic Forcep Curved, 1 x 2 Teeth

| Catalog # | Size         |
|-----------|--------------|
| MM-880-07 | 12 cm (4 ¾") |
| MM-880-08 | 14 cm (5 ½") |



### Halsted-Mosquito Heamostatic Forcep Straight

| Catalog # | Size           |
|-----------|----------------|
| MM-880-09 | 18.5 cm (7 ¼") |
| MM-880-10 | 21 cm (8 ¼")   |






### Halsted-Mosquito Hemostatic Forcep Curved

| Catalog # | Size         |
|-----------|--------------|
| MM-880-11 | 18 cm (7")   |
| MM-880-12 | 20.5 cm (8") |




### DeBakey-Mosquito Atrauma Hemostatic Forcep Straight

| Catalog # | Size       |
|-----------|------------|
| MM-880-13 | 15 cm (6") |



### DeBakey-Mosquito Atrauma Hemostatic Forcep Curved

| Catalog # | Size       |
|-----------|------------|
| MM-880-14 | 15 cm (6") |



### Crile Hemostatic Forcep Straight, 1 x 2 Teeth

| Catalog # | Size           |
|-----------|----------------|
| MM-880-15 | 14 cm (5 1/2") |



### Crile Hemostatic Forcep Curved, 1 x 2 Teeth

| Catalog # | Size           |
|-----------|----------------|
| MM-880-16 | 14 cm (5 1/2") |





### Spencer-Wells Haemostatic Forcep Straight

| Catalog # | Size           |
|-----------|----------------|
| MM-880-17 | 13 cm (5")     |
| MM-880-18 | 14.5 cm (5 ¾") |
| MM-880-19 | 15.5 cm (6")   |
| MM-880-20 | 18.5 cm (7 ¼") |
| MM-880-21 | 20 cm (8")     |




### Spencer-Wells Haemostatic Forcep Curved

| Catalog # | Size         |
|-----------|--------------|
| MM-880-22 | 13 cm (5")   |
| MM-880-23 | 14 cm (5 ¾") |
| MM-880-24 | 15 cm (6")   |
| MM-880-25 | 18 cm (7 ¼") |
| MM-880-26 | 20 cm (8")   |



### Pean Haemostatic Forcep

| Catalog # | Size           |
|-----------|----------------|
| MM-880-27 | 14.5 cm (5 ¾") |
| MM-880-28 | 16 cm (6 ¼")   |



### Allis-Baby Intestinal and Tissue Grasping Forcep 4 x 5 Teeth

| Catalog # | Size       |
|-----------|------------|
| MM-880-29 | 13 cm (5") |



### Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth

| Catalog # | Size         |
|-----------|--------------|
| MM-880-30 | 15.5 cm (6") |





### Young Tongue Holding Forcep

| Catalog # | Size       |
|-----------|------------|
| MM-880-31 | 15 cm (6") |



### Kelly Hemostatic Forcep Straight

| Catalog # | Size             |
|-----------|------------------|
| MM-880-32 | 14.5 cm (5 3/4") |



### Jones Towel Clamp Sharp

| Catalog # | Size            |
|-----------|-----------------|
| MM-880-33 | 5.5 cm (2 1/2") |
| MM-880-34 | 8.5 cm (3 1/4") |



### Lorna Towel Clamp Fine Serrated, Wide Jaw

| Catalog # | Size            |
|-----------|-----------------|
| MM-880-35 | 9.5 cm (3 3/4") |
| MM-880-36 | 13 cm (5")      |





✉ [info@makromed.tech](mailto:info@makromed.tech) | [sales@makromed.tech](mailto:sales@makromed.tech)

🌐 [www.makromed.tech](http://www.makromed.tech)

